

KD-Validated Anti-HRAS Rabbit Monoclonal Antibody

Rabbit monoclonal antibody


ABG1548

Product Overview


Name	KD-Validated Anti-HRAS Rabbit Monoclonal Antibody
Accession(Primary)	P01112
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human GTPase HRAS
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 21 kDa , observed , 19 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000 Flow Cytometry (FCM): 1:2,000

Target information (P01112)


Gene Name	HRAS
About Gene Symbol	HRAS (HRas Proto-Oncogene, GTPase) is a formin actin nucleator involved in cell signaling and proliferation. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Copyright ©2024 Genuin Biotechnologies LLC



Copyright ©2024 Genuin Biotechnologies LLC



Copyright ©2024 Genuin Biotechnologies LLC