

KD-Validated Anti-F11 receptor Rabbit Monoclonal Antibody

Rabbit monoclonal antibody


ABG1549

Product Overview


Name	KD-Validated Anti-F11 receptor Rabbit Monoclonal Antibody
Accession(Primary)	Q9Y624
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human JAM1
Species Reactivity	Human
Calculated MW	Predicted, 33 kDa, observed, 35 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000 Flow Cytometry (FCM): 1:2,000 Immunocytochemistry (IC): 1:1,000

Target information (Q9Y624)


Gene Name	F11R
About Gene Symbol	F11R (Junction adhesion molecule 1) is a membrane receptor involved in receptor-mediated signaling. It is an important research target in cell signaling, drug discovery, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.




Copyright ©2024 Genuin Biotechnologies LLC



Copyright ©2024 Genuin Biotechnologies LLC



Copyright ©2024 Genuin Biotechnologies LLC



Immunocytochemical staining of HAP-1 cells with F11 receptor antibody (Cat#ABG1549, 1:1,000). Nuclei were stained blue with DAPI; F11 receptor was stained magenta with Alexa Fluor® 647. Images were taken using Leica stellaris 5. Protein abundance based on laser Intensity and smart gain: Medium. Scale bar: 20 μ m.