

KD-Validated Anti-Keratin 10 Rabbit Monoclonal Antibody

Rabbit monoclonal antibody

ABG1550

Product Overview

Name	KD-Validated Anti-Keratin 10 Rabbit Monoclonal Antibody
Accession(Primary)	P13645
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human Cytokeratin 10
Species Reactivity	Human, Rat
Calculated MW	Predicted, 59 kDa , observed , 60 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000 Immunocytochemistry (IC): 1:100-1:1,000

Target information (P13645)

Gene Name	KRT10
About Gene Symbol	KRT10 (Keratin 10) is a keratin cytoskeletal protein involved in MAPK/ERK signaling. It is an important research target in cell biology, cancer metastasis, and cytoskeletal research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.

