

KD-Validated Anti-Histone deacetylase 6 Rabbit Monoclonal Antibody

Rabbit monoclonal antibody


ABG1552

Product Overview


Name	KD-Validated Anti-Histone deacetylase 6 Rabbit Monoclonal Antibody
Accession(Primary)	Q9UBN7
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human HDAC6
Species Reactivity	Human
Calculated MW	Predicted, 131 kDa , observed , 160 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:4,000-1:20,000 Flow Cytometry (FCM): 1:2,000

Target information (Q9UBN7)


Gene Name	HDAC6
About Gene Symbol	HDAC6 (Histone Deacetylase 6) is a histone deacetylase involved in chromatin remodeling and epigenetics. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Copyright ©2024 Genuin Biotechnologies LLC



Copyright ©2024 Genuin Biotechnologies LLC



Copyright ©2024 Genuin Biotechnologies LLC