

KD-Validated Anti-JUNB Rabbit Monoclonal Antibody

Rabbit monoclonal antibody




ABG1553


Product Overview

Name	KD-Validated Anti-JUNB Rabbit Monoclonal Antibody
Accession(Primary)	P17275
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human JunB
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 36 kDa , observed , 43 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000 Flow Cytometry (FCM): 1:2,000 Immunocytochemistry (IC): 1:1,000

Target information (P17275)

Gene Name	JUNB
About Gene Symbol	JUNB (JunB Proto-Oncogene, AP-1 Transcription Factor Subunit) is a proto-oncogene involved in gene transcription and regulation. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.





Immunocytochemical staining of HepG2 cells with anti-JUNB antibody (Cat#ABG1553, 1:1,000).

Nuclei were stained blue with DAPI; JUNB was stained magenta with Alexa Fluor® 647. Images were taken using Leica stellaris 5. Protein abundance based on laser Intensity and smart gain: Medium. Scale bar: 20 μ m.