

KD-Validated Anti-Hsp105 Rabbit Monoclonal Antibody

Rabbit monoclonal antibody

ABG1554

Product Overview

Name	KD-Validated Anti-Hsp105 Rabbit Monoclonal Antibody
Accession(Primary)	Q92598
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human Hsp105
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 97 kDa , observed , 105 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000 Flow Cytometry (FCM): 1:2,000

Target information (Q92598)

Gene Name	HSPH1
About Gene Symbol	HSPH1 (Heat Shock Protein Family H (Hsp110) Member 1) is a molecular chaperone involved in adaptive immunity. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.

