

KD-Validated Anti-Galectin 1 Rabbit Monoclonal Antibody

Rabbit monoclonal antibody


ABG1559

Product Overview


Name	KD-Validated Anti-Galectin 1 Rabbit Monoclonal Antibody
Accession(Primary)	P09382
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human Galectin 1
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 15 kDa , observed , 15 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000 Flow Cytometry (FCM): 1:2,000

Target information (P09382)


Gene Name	LGALS1
About Gene Symbol	LGALS1 (Galectin 1) is an extracellular matrix laminin protein involved in MAPK/ERK signaling. It is an important research target in cell biology, cancer metastasis, and cytoskeletal research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Copyright ©2024 Genuin Biotechnologies LLC



Copyright ©2024 Genuin Biotechnologies LLC



Copyright ©2024 Genuin Biotechnologies LLC