

KD-Validated Anti-Caspase 6 Rabbit Monoclonal Antibody

Rabbit monoclonal antibody


ABG1563

Product Overview


Name	KD-Validated Anti-Caspase 6 Rabbit Monoclonal Antibody
Accession(Primary)	P55212
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human Caspase-6
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 33 kDa , observed, 34 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB):1:1,000-1:5,000 Flow Cytometry (FCM): 1:2,000 Immunocytochemistry (IC): 1:100-1:1,000

Target information (P55212)


Gene Name	CASP6
About Gene Symbol	CASP6 (Caspase 6) is a cysteine protease (caspase) involved in apoptosis. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.




Western blotting analysis using anti-Caspase 6 antibody (Cat#ABG1563). Total cell lysates (30 µg) from various cell lines were loaded and separated by SDS-PAGE. The blot was incubated with anti-Caspase 6 antibody (Cat#ABG1563, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody respectively.



Western blotting analysis using anti-caspase 6 antibody (Cat#ABG1563). Caspase 6 expression in wild-type (WT) and caspase 6 (CASP6) shRNA knockdown (KD) C2C12 cells with 20 µg of total cell lysates. β-Tubulin serves as a loading control. The blot was incubated with anti-caspase 6 antibody (Cat#ABG1563, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody respectively.



Flow cytometric analysis of Caspase 6 expression in HT-1080 cells using Caspase 6 antibody (Cat#ABG1563, 1:2,000). Green, isotype control; red, Caspase 6.



Immunocytochemical staining of HT-1080 cells with Caspase 6 antibody (Cat#ABG1563, 1:1,000). Nuclei were stained blue with DAPI; Caspase 6 was stained magenta with Alexa Fluor® 647. Images were taken using Leica stellaris 5. Protein abundance based on laser Intensity and smart gain: Medium. Scale bar: 20 μ m.