

KD-Validated Anti-Goosecoid homeobox Rabbit Monoclonal Antibody

Rabbit monoclonal antibody


ABG1568

Product Overview

Name	KD-Validated Anti-Goosecoid homeobox Rabbit Monoclonal Antibody
Accession(Primary)	P56915
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human GSC
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 28 kDa , observed , 28 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB):1:1,000-1:5,000


Target information (P56915)

Gene Name	GSC
About Gene Symbol	GSC (Goosecoid Homeobox) is a transcription factor involved in gene transcription and regulation. It is an important research target in gene regulation, development, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Copyright ©2024 Genuin Biotechnologies LLC

Copyright ©2024 Genuin Biotechnologies LLC



Western blotting analysis using anti-Goosecoid homeobox antibody (Cat#ABG1568). Total cell lysates (30 µg) from various cell lines were loaded and separated by SDS-PAGE. The blot was incubated with anti-Goosecoid homeobox antibody (Cat#ABG1568, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody respectively.

Western blotting analysis using anti-Goosecoid homeobox antibody (Cat#ABG1568). Goosecoid homeobox expression in wild type (WT) and Goosecoid homeobox shRNA knockdown (KD) HeLa cells with 30 µg of total cell lysates. β-Tubulin serves as a loading control. The blot was incubated with anti-Goosecoid homeobox antibody (Cat#ABG1568, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody respectively.