

KD-Validated Anti-Phospho-JunD (S255) Rabbit Monoclonal Antibody

Rabbit monoclonal antibody


ABG1570

Product Overview


Name	KD-Validated Anti-Phospho-JunD (S255) Rabbit Monoclonal Antibody
Accession(Primary)	P17535
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human JunD
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 35 kDa , observed , 38,42 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB):1:1,000-1:5,000 Flow Cytometry (FCM): 1:2,000 Immunocytochemistry (IC): 1:1,000

Target information (P17535)


Gene Name	JUND
About Gene Symbol	JUND (JunD Proto-Oncogene, AP-1 Transcription Factor Subunit) is a proto-oncogene involved in gene transcription and regulation. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.




Copyright ©2024 Genuin Biotechnologies LLC



Copyright ©2024 Genuin Biotechnologies LLC



Copyright ©2024 Genuin Biotechnologies LLC



Immunocytochemical staining of HepG2 cells with anti-Phospho-JunD (S255) antibody (Cat#ABG1570, 1:1,000). Nuclei were stained blue with DAPI; Phospho-JunD (S255) was stained magenta with Alexa Fluor® 647. Images were taken using Leica stellaris 5. Protein abundance based on laser Intensity and smart gain: Medium. Scale bar: 20 μ m.