

KD-Validated Anti-Glutamic-oxaloacetic transaminase 2 Rabbit Monoclonal Antibody

Rabbit monoclonal antibody


ABG1579

Product Overview


Name	KD-Validated Anti-Glutamic-oxaloacetic transaminase 2 Rabbit Monoclonal Antibody
Accession(Primary)	P00505
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human GOT2
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 48 kDa , observed , 43 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB):1:1,000-1:5,000 Flow Cytometry (FCM):1:2,000

Target information (P00505)


Gene Name	GOT2
About Gene Symbol	GOT2 (Glutamic-Oxaloacetic Transaminase 2) is a membrane receptor involved in receptor-mediated signaling. It is an important research target in cell signaling, drug discovery, and cancer research and is widely studied to elucidate cellular mechanisms and



Copyright ©2024, Genuin Biotechnologies LLC



Copyright ©2024, Genuin Biotechnologies LLC



Copyright ©2024, Genuin Biotechnologies LLC