

KD-Validated Anti-Mevalonate Kinase Rabbit Monoclonal Antibody

Rabbit monoclonal antibody

ABG1588

Product Overview

Name	KD-Validated Anti-Mevalonate Kinase Rabbit Monoclonal Antibody
Accession(Primary)	Q03426
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human MVK
Species Reactivity	Human, Rat
Calculated MW	Predicted, 42 kDa, observed, 42 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB):1:1,000-1:5,000 Immunocytochemistry (IC): 1:100-1:1,000

Target information (Q03426)

Gene Name	MVK
About Gene Symbol	MVK (Mevalonate Kinase) is a membrane receptor involved in translational regulation. It is an important research target in cell signaling, drug discovery, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.

