

KD-Validated Anti-Major Vault Protein Rabbit Monoclonal Antibody

Rabbit monoclonal antibody


ABG1607

Product Overview


Name	KD-Validated Anti-Major Vault Protein Rabbit Monoclonal Antibody
Accession(Primary)	Q14764
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human MVP
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 99 kDa , observed , 110 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB):1:1,000-1:5,000 Flow Cytometry (FCM): 1:2,000 Immunocytochemistry (IC): 1:100-1:1,000

Target information (Q14764)


Gene Name	MVP
About Gene Symbol	MVP (Major Vault Protein) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.




Copyright ©2024 Genuin Biotechnologies LLC



Copyright ©2024 Genuin Biotechnologies LLC



Copyright ©2024 Genuin Biotechnologies LLC



Immunocytochemical staining of HeLa cells with anti-Major vault protein antibody (Cat#ABG1607, 1:1,000). Nuclei were stained blue with DAPI; Major vault protein was stained magenta with Alexa Fluor® 647. Images were taken using Leica stellaris 5. Protein abundance based on laser Intensity and smart gain: Medium. Scale bar: 20 μ m.