

KD-Validated Anti-NSDHL Rabbit Monoclonal Antibody

Rabbit monoclonal antibody




ABG1610


Product Overview

Name	KD-Validated Anti-NSDHL Rabbit Monoclonal Antibody
Accession(Primary)	Q15738
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human NSDHL
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 42 kDa , observed , 38 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB):1:1,000-1:5,000 Immunocytochemistry (IC): 1:100-1:1,000 Flow Cytometry (FCM):1:2,000

Target information (Q15738)

Gene Name	NSDHL
About Gene Symbol	NSDHL (NAD(P) Dependent Steroid Dehydrogenase-Like) is a metabolic dehydrogenase involved in cellular metabolism. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular me





Immunocytochemical staining of HT-1080 cells with anti-NSDHL antibody (Cat#ABG1610, 1:1,000). Nuclei were stained blue with DAPI; NSDHL was stained magenta with Alexa Fluor® 647. Images were taken using Leica stellaris 5. Protein abundance based on laser Intensity and smart gain: Medium. Scale bar: 20 μ m.