

KD-Validated Anti-Nemo Like Kinase Rabbit Monoclonal Antibody

Rabbit monoclonal antibody


ABG1616

Product Overview


Name	KD-Validated Anti-Nemo Like Kinase Rabbit Monoclonal Antibody
Accession(Primary)	Q9UBE8
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human NLK
Species Reactivity	Human
Calculated MW	Predicted, 58 kDa ; Observed , 58 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB):1:1,000-1:5,000 Flow Cytometry (FCM): 1:2,000 Immunocytochemistry (IC): 1:100-1:1,000

Target information (Q9UBE8)


Gene Name	NLK
About Gene Symbol	NLK (Nemo Like Kinase) is a serine/threonine protein kinase involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.




Copyright ©2024 Genuin Biotechnologies LLC



Copyright ©2024 Genuin Biotechnologies LLC



Copyright ©2024 Genuin Biotechnologies LLC



Immunocytochemical staining of HeLa cells with anti-Nemo like kinase antibody (Cat#ABG1616, 1:1,000). Nuclei were stained blue with DAPI; Nemo like kinase was stained magenta with Alexa Fluor® 647. Images were taken using Leica stellaris 5. Protein abundance based on laser Intensity and smart gain: Medium. Scale bar: 20 μ m.