

KD-Validated Anti-NDUFA13 Rabbit Monoclonal Antibody

Rabbit monoclonal antibody

ABG1617

Product Overview

Name	KD-Validated Anti-NDUFA13 Rabbit Monoclonal Antibody
Accession(Primary)	Q9P0J0
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human GRIM19
Species Reactivity	Human
Calculated MW	Predicted, 17 kDa , observed , 17 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB):1:1,000-1:5,000 Immunocytochemistry (IC): 1:100-1:1,000

Target information (Q9P0J0)

Gene Name	NDUFA13
About Gene Symbol	NDUFA13 (NADH:Ubiquinone Oxidoreductase Subunit A13) is a metabolic dehydrogenase involved in JAK-STAT/cytokine signaling and innate immune signaling. It is an important research target in immunology, inflammation, and cancer research and is widely studie

