

KD-Validated Anti-LYN Proto-Oncogene Rabbit Monoclonal Antibody

Rabbit monoclonal antibody




ABG1621


Product Overview

Name	KD-Validated Anti-LYN Proto-Oncogene Rabbit Monoclonal Antibody
Accession(Primary)	P07948
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human Lyn
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 59 kDa , observed , 56 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000 Immunocytochemistry (IC): 1:100-1:1,000 Flow Cytometry (FCM):1:2,000

Target information (P07948)

Gene Name	LYN
About Gene Symbol	LYN (LYN Proto-Oncogene, Src Family Tyrosine Kinase) is a receptor tyrosine kinase involved in cell signaling and proliferation. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.





Immunocytochemical staining of HepG2 cells with anti-LYN antibody (Cat#ABG1621, 1:1,000). Nuclei were stained blue with DAPI; LYN was stained magenta with Alexa Fluor® 647. Images were taken using Leica stellaris 5. Protein abundance based on laser Intensity and smart gain: High. Scale bar: 20 μ m.