

KD-Validated Anti-Pyruvate Dehydrogenase E2 Rabbit Monoclonal Antibody

Rabbit monoclonal antibody


ABG1623

Product Overview


Name	KD-Validated Anti-Pyruvate Dehydrogenase E2 Rabbit Monoclonal Antibody
Accession(Primary)	P10515
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human Pyruvate Dehydrogenase E2
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted , 69 kDa , observed , 69 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000 Flow Cytometry (FCM): 1:2,000 Immunocytochemistry (IC): 1:1,000

Target information (P10515)


Gene Name	DLAT
About Gene Symbol	DLAT (Dihydrolipoamide S-Acetyltransferase) is a metabolic dehydrogenase involved in adaptive immunity. It is an important research target in fibrotic disease and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.




Copyright ©2024 Genuin Biotechnologies LLC



Copyright ©2024 Genuin Biotechnologies LLC



Copyright ©2024 Genuin Biotechnologies LLC



Immunocytochemical staining of HepG2 cells with Pyruvate Dehydrogenase E2 antibody (Cat#ABG1623, 1:1,000). Nuclei were stained blue with DAPI; Pyruvate Dehydrogenase E2 was stained magenta with Alexa Fluor® 647. Images were taken using Leica stellaris 5. Protein abundance based on laser Intensity and smart gain: Low. Scale bar: 20 μ m.