

KD-Validated Anti-Phospho-Rad17 (S656) Rabbit Monoclonal Antibody

Rabbit monoclonal antibody


ABG1626

Product Overview


Name	KD-Validated Anti-Phospho-Rad17 (S656) Rabbit Monoclonal Antibody
Accession(Primary)	O75943
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human Phospho-Rad17 (S656)
Species Reactivity	Human, Mouse
Calculated MW	Predicted, 77 kDa , observed , 85 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB):1:1,000-1:5,000 Flow Cytometry(FCM):1:2,000

Target information (O75943)


Gene Name	RAD17
About Gene Symbol	RAD17 (RAD17 Checkpoint Clamp Loader Component) is a protein kinase and signaling regulator involved in cell cycle regulation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Copyright ©2024 Genuin Biotechnologies LLC



Copyright ©2024 Genuin Biotechnologies LLC



Copyright ©2024 Genuin Biotechnologies LLC