

KD-Validated Anti-RAVER2 Rabbit Monoclonal Antibody

Rabbit monoclonal antibody




ABG1627


Product Overview

Name	KD-Validated Anti-RAVER2 Rabbit Monoclonal Antibody
Accession(Primary)	Q9HCJ3
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human RAVER2
Species Reactivity	Human
Calculated MW	Predicted, 74 kDa , observed , 74 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB):1:1,000-1:5,000 Immunocytochemistry (IC): 1:100-1:1,000 Flow Cytometry (FCM): 1:2,000

Target information (Q9HCJ3)

Gene Name	RAVER2
About Gene Symbol	RAVER2 (Ribonucleoprotein, PTB Binding 2) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.





Immunocytochemical staining of HeLa cells with anti-RAVER2 antibody (Cat#ABG1627, 1:1,000). Nuclei were stained blue with DAPI; RAVER2 was stained magenta with Alexa Fluor® 647. Images were taken using Leica stellaris 5. Protein abundance based on laser Intensity and smart gain: Medium. Scale bar: 20 μ m.