

KD-Validated Anti-RelB Rabbit Monoclonal Antibody

Rabbit monoclonal antibody




ABG1637


Product Overview

Name	KD-Validated Anti-RelB Rabbit Monoclonal Antibody
Accession(Primary)	Q01201
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human RelB
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 62 kDa , observed , 70 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB):1:1,000-1:5,000 Immunocytochemistry (IC): 1:100-1:1,000 Flow Cytometry (FCM):1:2,000

Target information (Q01201)

Gene Name	RELB
About Gene Symbol	RELB (RELB Proto-Oncogene, NF- κ B Subunit) is a proto-oncogene involved in NF- κ B signaling and angiogenesis. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.





Immunocytochemical staining of C2C12 cells with anti-RELB antibody (Cat#ABG1637, 1:1,000).

Nuclei were stained blue with DAPI; RELB was stained magenta with Alexa Fluor® 647. Images were taken using Leica stellaris 5. Protein abundance based on laser Intensity and smart gain: Medium. Scale bar: 20 μ m.