

KD-Validated Anti-Retinoic Acid Receptor beta Rabbit Monoclonal Antibody

Rabbit monoclonal antibody


ABG1647

Product Overview


Name	KD-Validated Anti-Retinoic Acid Receptor beta Rabbit Monoclonal Antibody
Accession(Primary)	P10826
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human Retinoic Acid Receptor beta
Species Reactivity	Human
Calculated MW	Predicted, 50 kDa ; Observed , 55 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB):1:1,000-1:5,000 Immunocytochemistry (IC): 1:100-1:1,000 Flow Cytometry (FCM): 1:2,000

Target information (P10826)


Gene Name	RARB
About Gene Symbol	RARB (Retinoic Acid Receptor Beta) is a membrane receptor involved in receptor-mediated signaling. It is an important research target in infectious disease and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.




Copyright ©2024 Genuin Biotechnologies LLC



Copyright ©2024 Genuin Biotechnologies LLC



Copyright ©2024 Genuin Biotechnologies LLC



Immunocytochemical staining of HeLa cells with anti-Retinoic acid receptor beta antibody (Cat#ABG1647, 1:1,000). Nuclei were stained blue with DAPI; Retinoic acid receptor beta was stained magenta with Alexa Fluor® 647. Images were taken using Leica stellaris 5. Protein abundance based on laser Intensity and smart gain: Medium. Scale bar: 20 μ m.