

KD-Validated Anti-Ribophorin I Rabbit Monoclonal Antibody

Rabbit monoclonal antibody

ABG1651

Product Overview

Name	KD-Validated Anti-Ribophorin I Rabbit Monoclonal Antibody
Accession(Primary)	P04843
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human PRibophorin I
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 69 kDa , observed , 69 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB):1:1,000-1:5,000 Flow Cytometry (FCM):1:2,000

Target information (P04843)

Gene Name	RPN1
About Gene Symbol	RPN1 (Ribophorin I) is a metabolic enzyme involved in cellular metabolism. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.

