

KD-Validated Anti-Ring Finger Protein 40 Rabbit Monoclonal Antibody

Rabbit monoclonal antibody

ABG1653

Product Overview

Name	KD-Validated Anti-Ring Finger Protein 40 Rabbit Monoclonal Antibody
Accession(Primary)	O75150
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human RNF40
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 114 kDa , observed , 130 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB):1:1,000-1:5,000 Immunocytochemistry (IC): 1:100-1:1,000

Target information (O75150)

Gene Name	RNF40
About Gene Symbol	RNF40 (Ring Finger Protein 40) is a ubiquitin E3 ligase involved in ubiquitin-proteasome pathway. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.

