

KD-Validated Anti-PELP1 Rabbit Monoclonal Antibody

Rabbit monoclonal antibody

ABG1655

Product Overview

Name	KD-Validated Anti-PELP1 Rabbit Monoclonal Antibody
Accession(Primary)	Q8IZL8
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human PELP1
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 120 kDa , observed , 160 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB):1:1,000-1:5,000 Flow Cytometry (FCM): 1:2,000

Target information (Q8IZL8)

Gene Name	PELP1
About Gene Symbol	PELP1 (Proline, Glutamate And Leucine Rich Protein 1) is a membrane receptor involved in receptor-mediated signaling. It is an important research target in cell signaling, drug discovery, and cancer research and is widely studied to elucidate cellular mec

