

KD-Validated Anti-Retinoic Acid Receptor Alpha Rabbit

Rabbit monoclonal antibody


ABG1657

Product Overview


Name	KD-Validated Anti-Retinoic Acid Receptor Alpha Rabbit Monoclonal Antibody
Accession(Primary)	P10276
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human Retinoic Acid Receptor alpha
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 51 kDa , o bserved , 55 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB):1:1,000-1:5,000 Immunocytochemistry (IC): 1:100-1:1,000

Target information (P10276)


Gene Name	RARA
About Gene Symbol	RARA (Retinoic Acid Receptor Alpha) is a membrane receptor involved in receptor-mediated signaling. It is an important research target in cell signaling, drug discovery, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Copyright ©2024 Genuin Biotechnologies LLC



Copyright ©2024 Genuin Biotechnologies LLC



Copyright ©2024 Genuin Biotechnologies LLC