

KD-Validated Anti-Phospho-PLCB3 (S537) Rabbit Monoclonal Antibody

Rabbit monoclonal antibody

ABG1667

Product Overview

Name	KD-Validated Anti-Phospho-PLCB3 (S537) Rabbit Monoclonal Antibody
Accession(Primary)	Q01970
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human Phospho-PLCB3 (S537)
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 139 kDa , observed , 150 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB):1:1,000-1:5,000

Target information (Q01970)

Gene Name	PLCB3
About Gene Symbol	PLCB3 (Phospholipase C Beta 3) is a metabolic enzyme involved in cellular metabolism. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.

