

KD-Validated Anti-PRKAR2A Rabbit Monoclonal Antibody

Rabbit monoclonal antibody


ABG1668

Product Overview


Name	KD-Validated Anti-PRKAR2A Rabbit Monoclonal Antibody
Accession(Primary)	P13861
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human PKA R2
Species Reactivity	Human
Calculated MW	Predicted, 46 kDa , observed , 52 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB):1:1,000-1:5,000 Flow Cytometry(FCM):1:2,000

Target information (P13861)


Gene Name	PRKAR2A
About Gene Symbol	PRKAR2A (Protein Kinase CAMP-Dependent Type II Regulatory Subunit Alpha) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research



Copyright ©2024 Genuin Biotechnologies LLC



Copyright ©2024 Genuin Biotechnologies LLC



Copyright ©2024 Genuin Biotechnologies LLC