

KD-Validated Anti-NHERF1 Rabbit Monoclonal Antibody

Rabbit monoclonal antibody


ABG1669

Product Overview


Name	KD-Validated Anti-NHERF1 Rabbit Monoclonal Antibody
Accession(Primary)	O14745
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human EBP50
Species Reactivity	Human, Rat
Calculated MW	Predicted, 39 kDa , observed , 50 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB):1:1,000-1:5,000 Flow Cytometry (FCM): 1:2,000

Target information (O14745)


Gene Name	NHERF1
About Gene Symbol	NHERF1 (NHERF Family PDZ Scaffold Protein 1) is a signaling scaffold protein involved in membrane transport. It is an important research target in pharmacology, drug resistance, and membrane biology and is widely studied to elucidate cellular mechanisms a



Copyright ©2024 Genuin Biotechnologies LLC



Copyright ©2024 Genuin Biotechnologies LLC



Copyright ©2024 Genuin Biotechnologies LLC

