

KD-Validated Anti-Phospholipid Scramblase 3 Rabbit Monoclonal Antibody

Rabbit monoclonal antibody


ABG1675

Product Overview


Name	KD-Validated Anti-Phospholipid Scramblase 3 Rabbit Monoclonal Antibody
Accession(Primary)	Q9NRY6
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human PLSCR3
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 32 kDa , observed, 32 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB):1:1,000-1:5,000 Flow Cytometry (FCM):1:2,000

Target information (Q9NRY6)


Gene Name	PLSCR3
About Gene Symbol	PLSCR3 (Phospholipid Scramblase 3) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Copyright ©2024 Genuin Biotechnologies LLC



right ©2024 Genuin Biotechnologies L



Copyright ©2024 Genuin Biotechnologies LLC