

KD-Validated Anti-Inorganic Pyrophosphatase 1 Rabbit

Rabbit monoclonal antibody


ABG1676

Product Overview


Name	KD-Validated Anti-Inorganic Pyrophosphatase 1 Rabbit Monoclonal Antibody
Accession(Primary)	Q15181
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human Pyrophosphatase 1
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 33 kDa , observed, 33 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB):1:1,000-1:5,000 Flow Cytometry (FCM):1:2,000

Target information (Q15181)


Gene Name	PPA1
About Gene Symbol	PPA1 (Inorganic Pyrophosphatase 1) is a protein phosphatase involved in cellular metabolism. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Copyright ©2024 Genuin Biotechnologies LLC



Copyright ©2024 Genuin Biotechnologies LLC



Copyright ©2024 Genuin Biotechnologies LLC