

KD-Validated Anti-Pancreatic And Duodenal Homeobox 1 Rabbit Monoclonal Antibody

Rabbit monoclonal antibody

ABG1677

Product Overview

Name	KD-Validated Anti-Pancreatic And Duodenal Homeobox 1 Rabbit Monoclonal Antibody
Accession(Primary)	P52945
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human PDX1
Species Reactivity	Human, Rat
Calculated MW	Predicted, 31 kDa , observed, 42 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB):1:1,000-1:5,000 Flow Cytometry (FCM): 1:2,000

Target information (P52945)

Gene Name	PDX1
About Gene Symbol	PDX1 (Pancreatic And Duodenal Homeobox 1) is a transcription factor involved in gene transcription and regulation. It is an important research target in gene regulation, development, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



