

KD-Validated Anti-Prostaglandin E Synthase 3 Rabbit Monoclonal Antibody


ABG1681

Product Overview


Name	KD-Validated Anti-Prostaglandin E Synthase 3 Rabbit Monoclonal Antibody
Accession(Primary)	Q15185
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human p23
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 19 kDa , observed , 23 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western blotting (WB):1:1,000-1:5,000 Flow Cytometry (FCM):1:2,000

Target information (Q15185)


Gene Name	PTGES3
About Gene Symbol	PTGES3 (Prostaglandin E Synthase 3) is a biosynthetic enzyme involved in receptor-mediated signaling. It is an important research target in cell signaling, drug discovery, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Copyright ©2024 Genuin Biotechnologies LLC



Copyright ©2024 Genuin Biotechnologies LLC



Copyright ©2024 Genuin Biotechnologies LLC