

KD-Validated Anti-SCARB2 Rabbit Monoclonal Antibody

Rabbit monoclonal antibody


ABG1687

Product Overview


Name	KD-Validated Anti-SCARB2 Rabbit Monoclonal Antibody
Accession(Primary)	Q14108
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human LIMPII
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 54 kDa , observed , 77 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000 Flow Cytometry (FCM):1:2,000

Target information (Q14108)


Gene Name	SCARB2
About Gene Symbol	SCARB2 (Scavenger Receptor Class B Member 2) is a collagen extracellular matrix protein involved in adaptive immunity and vesicular trafficking. It is an important research target in cell signaling, drug discovery, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Copyright ©2024 Genuin Biotechnologies LLC



Copyright ©2024 Genuin Biotechnologies LLC



Copyright ©2024 Genuin Biotechnologies LLC