

KD-Validated Anti-PIWIL1 Rabbit Monoclonal Antibody

Rabbit monoclonal antibody

ABG1690

Product Overview

Name	KD-Validated Anti-PIWIL1 Rabbit Monoclonal Antibody
Accession(Primary)	Q96J94
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human PIWIL1
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 99 kDa , observed, 99 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB):1:1,000-1:5,000 Flow Cytometry (FCM):1:2,000

Target information (Q96J94)

Gene Name	PIWIL1
About Gene Symbol	PIWIL1 (Piwi Like RNA-Mediated Gene Silencing 1) is a metabolic enzyme involved in cellular metabolism. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.

