

KD-Validated Anti-Phospholipid Scramblase 1 Rabbit Monoclonal Antibody

Rabbit monoclonal antibody


ABG1691

Product Overview


Name	KD-Validated Anti-Phospholipid Scramblase 1 Rabbit Monoclonal Antibody
Accession(Primary)	O15162
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human Scramblase 1
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 35 kDa , observed, 35 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB):1:1,000-1:5,000 Flow Cytometry (FCM):1:2,000

Target information (O15162)


Gene Name	PLSCR1
About Gene Symbol	PLSCR1 (Phospholipid Scramblase 1) is a metabolic enzyme involved in cellular metabolism. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Copyright © 2024 Genuin Biotechnologies LLC



Copyright © 2024 Genuin Biotechnologies LLC



Copyright © 2024 Genuin Biotechnologies LLC