

KD-Validated Anti-PERP Rabbit Monoclonal Antibody

Rabbit monoclonal antibody


ABG1697

Product Overview


Name	KD-Validated Anti-PERP Rabbit Monoclonal Antibody
Accession(Primary)	Q96FX8
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human PERP
Species Reactivity	Human
Calculated MW	Predicted, 21 kDa , observed , 21 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000 Flow Cytometry (FCM):1:2,000

Target information (Q96FX8)


Gene Name	PERP
About Gene Symbol	PERP (P53 Apoptosis Effector Related To PMP22) is a keratin cytoskeletal protein involved in apoptosis. It is an important research target in cell signaling, drug discovery, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Copyright © 2024 Genuin Biotechnologies LLC



ht © 2024 Genuin Biotechnologies LLC



Copyright © 2024 Genuin Biotechnologies LLC