

KD-Validated Anti-Phosphoglycerate Mutase 1 Rabbit

Rabbit monoclonal antibody


ABG1698

Product Overview


Name	KD-Validated Anti-Phosphoglycerate Mutase 1 Rabbit Monoclonal Antibody
Accession(Primary)	P18669
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human PGAM1
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 29 kDa , o bserved , 29 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB):1:1,000-1:5,000

Target information (P18669)

Gene Name	PGAM1
About Gene Symbol	PGAM1 (Phosphoglycerate Mutase 1) is a metabolic enzyme involved in cellular metabolism. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Copyright ©2024 Genulin Biotechnologies LLC



Copyright ©2024 Genulin Biotechnologies LLC