

KD-Validated Anti-Ubiquitin Specific Peptidase 39 Rabbit Monoclonal Antibody

Rabbit monoclonal antibody


ABG1700

Product Overview


Name	KD-Validated Anti-Ubiquitin Specific Peptidase 39 Rabbit Monoclonal Antibody
Accession(Primary)	Q53GS9
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human USP39
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 65 kDa, observed , 65 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000 Immunocytochemistry (IC): 1:100-1:1,000 Flow Cytometry (FCM): 1:2,000

Target information (Q53GS9)


Gene Name	USP39
About Gene Symbol	USP39 (Ubiquitin Specific Peptidase 39) is a protease enzyme involved in ubiquitin-proteasome pathway. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.




Copyright ©2025 Genuin Biotechnologies LLC



©2024 Genuin Biotechnologies LLC



Copyright ©2024 Genuin Biotechnologies LLC



Immunocytochemical staining of HepG2 cells with anti-USP39 antibody (Cat#ABG1700, 1:1,000).

Nuclei were stained blue with DAPI; USP39 was stained magenta with Alexa Fluor® 647. Images were taken using Leica stellaris 5. Protein abundance based on laser Intensity and smart gain: Medium. Scale bar: 20 μ m.