

KD-Validated Anti-Stomatin Like 2 Rabbit Monoclonal Antibody

Rabbit monoclonal antibody


ABG1703

Product Overview


Name	KD-Validated Anti-Stomatin Like 2 Rabbit Monoclonal Antibody
Accession(Primary)	Q9UJZ1
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human Stomatin like 2
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 39 kDa, observed, 36 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000 Immunocytochemistry (IC): 1:100-1:1,000 Flow Cytometry (FCM): 1:2,000

Target information (Q9UJZ1)


Gene Name	STOML2
About Gene Symbol	STOML2 (Stomatin Like 2) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.




Copyright ©2024 Genuin Biotechnologies LLC



©2024 Genuin Biotechnologies LLC



Copyright ©2024 Genuin Biotechnologies LLC



Immunocytochemical staining of HepG2 cells with anti-Stomatin like 2 antibody (Cat#ABG1703, 1:1,000). Nuclei were stained blue with DAPI; Stomatin like 2 was stained magenta with Alexa Fluor® 647. Images were taken using Leica stellaris 5. Protein abundance based on laser Intensity and smart gain: Medium. Scale bar: 20 μ m.