

KD-Validated Anti-Transcription Factor 7 Like 1 Rabbit Monoclonal Antibody

Rabbit monoclonal antibody


ABG1706

Product Overview


Name	KD-Validated Anti-Transcription Factor 7 Like 1 Rabbit Monoclonal Antibody
Accession(Primary)	Q9HCS4
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human TCF3
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 63 kDa ; Observed , 63 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000 Flow Cytometry (FCM): 1:2,000 Immunocytochemistry (IC): 1:100-1:1,000

Target information (Q9HCS4)


Gene Name	TCF7L1
About Gene Symbol	TCF7L1 (Transcription Factor 7 Like 1) is a transcription factor involved in gene transcription and regulation. It is an important research target in gene regulation, development, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.




Copyright ©2024 Genuin Biotechnologies LLC



©2024 Genuin Biotechnologies LLC



Copyright ©2024 Genuin Biotechnologies LLC



Immunocytochemical staining of HT-1080 cells with anti-transcription factor 7 like 1 antibody (Cat#ABG1706, 1:1,000). Nuclei were stained blue with DAPI; Transcription factor 7 like 1 was stained magenta with Alexa Fluor® 647. Images were taken using Leica stellaris 5. Protein abundance based on laser Intensity and smart gain: Medium. Scale bar: 20 μ m.