

KD-Validated Anti-TRK fused gene Rabbit Monoclonal Antibody

Rabbit monoclonal antibody




ABG1709


Product Overview

Name	KD-Validated Anti-TRK fused gene Rabbit Monoclonal Antibody
Accession(Primary)	Q92734
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human TRK fused gene
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 43 kDa , observed , 55 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:4,000-1:20,000 Immunocytochemistry (IC): 1:100-1:1,000 Flow Cytometry (FCM): 1:2,000

Target information (Q92734)

Gene Name	TFG
About Gene Symbol	TFG (Trafficking From ER To Golgi Regulator) is a protein kinase and signaling regulator involved in vesicular trafficking. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.





Immunocytochemical staining of HepG2 cells with anti-TFG antibody (Cat#ABG1709, 1:1,000). Nuclei were stained blue with DAPI; TFG was stained magenta with Alexa Fluor® 647. Images were taken using Leica stellaris 5. Protein abundance based on laser Intensity and smart gain: High. Scale bar: 20 μ m.