

KD-Validated Anti-UQCRC2 Rabbit Monoclonal Antibody

Rabbit monoclonal antibody




ABG1711


Product Overview

Name	KD-Validated Anti-UQCRC2 Rabbit Monoclonal Antibody
Accession(Primary)	P22695
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human UQCRC2
Species Reactivity	Human, Rat
Calculated MW	Predicted, 48 kDa, observed, 48 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB):1:1,000-1:5,000 Immunocytochemistry (IC): 1:100-1:1,000 Flow Cytometry (FCM): 1:2,000

Target information (P22695)

Gene Name	UQCRC2
About Gene Symbol	UQCRC2 (Ubiquinol-Cytochrome C Reductase Core Protein 2) is a metabolic enzyme involved in cellular metabolism. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.





Immunocytochemical staining of HepG2 cells with anti-UQCRC2 antibody (Cat#ABG1711, 1:1,000). Nuclei were stained blue with DAPI; UQCRC2 was stained magenta with Alexa Fluor® 647. Images were taken using Leica stellaris 5. Protein abundance based on laser Intensity and smart gain: Medium. Scale bar: 20 μ m.