

KD-Validated Anti-Tubulin Beta 3 Class III Rabbit Monoclonal Antibody

Rabbit monoclonal antibody

ABG1713

Product Overview

Name	KD-Validated Anti-Tubulin Beta 3 Class III Rabbit Monoclonal Antibody
Accession(Primary)	Q13509
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human beta Tubulin III
Species Reactivity	Human, Mouse
Calculated MW	Predicted, 50 kDa , observed , 55 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:2,000-1:10,000 Flow Cytometry (FCM): 1:2,000

Target information (Q13509)

Gene Name	TUBB3
About Gene Symbol	TUBB3 (Tubulin Beta 3 Class III) is a structural cytoskeletal protein involved in cytoskeletal organization. It is an important research target in fibrotic disease and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



