

KD-Validated Anti-VPS28 Subunit Of ESCRT-I Rabbit Monoclonal Antibody

Rabbit monoclonal antibody


ABG1715

Product Overview


Name	KD-Validated Anti-VPS28 Subunit Of ESCRT-I Rabbit Monoclonal Antibody
Accession(Primary)	Q9UK41
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human VPS28
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 25 kDa , observed , 25 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB):1:1,000-1:5,000 Flow Cytometry (FCM):1:2,000

Target information (Q9UK41)


Gene Name	VPS28
About Gene Symbol	VPS28 (VPS28 Subunit Of ESCRT-I) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cel



Western blotting analysis using anti-VPS28 antibody (Cat#ABG1715). Total cell lysates (30 μ g) from various cell lines were loaded and separated by SDS-PAGE. The blot was incubated with anti-VPS28 antibody (Cat#ABG1715, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody respectively.



Western blotting analysis using anti-VPS28 antibody (Cat#ABG1715). VPS28 expression in wild type (WT) and VPS28 shRNA knockdown (KD) HeLa cells with 20 μ g of total cell lysates. β -Tubulin serves as a loading control. The blot was incubated with anti-VPS28 antibody (Cat#ABG1715, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody respectively.



Flow cytometric analysis of VPS28 expression in C2C12 cells using anti-VPS28 antibody (Cat#ABG1715, 1:2,000). Green, isotype control; red, VPS28.

