

KD-Validated Anti-Ubiquitination Factor E4B Rabbit Monoclonal Antibody

Rabbit monoclonal antibody




ABG1716


Product Overview

Name	KD-Validated Anti-Ubiquitination Factor E4B Rabbit Monoclonal Antibody
Accession(Primary)	O95155
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human UBE4B
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 146 kDa , observed , 130 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB):1:1,000-1:5,000 Flow Cytometry (FCM): 1:2,000 Immunocytochemistry (IC): 1:100-1:1,000

Target information (O95155)

Gene Name	UBE4B
About Gene Symbol	UBE4B (Ubiquitination Factor E4B) is a metabolic enzyme involved in ubiquitin-proteasome pathway. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.





Immunocytochemical staining of HepG2 cells with anti-Ubiquitination factor E4 B antibody (Cat#ABG1716, 1:1,000). Nuclei were stained blue with DAPI; Ubiquitination factor E4 B was stained magenta with Alexa Fluor® 647. Images were taken using Leica stellaris 5. Protein abundance based on laser Intensity and smart gain: Medium. Scale bar: 20 μ m.