

KD-Validated Anti-Ubiquitination Factor E4A Rabbit Monoclonal Antibody


ABG1717

Product Overview


Name	KD-Validated Anti-Ubiquitination Factor E4A Rabbit Monoclonal Antibody
Accession(Primary)	Q14139
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human UBE4A
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 123 kDa , observed , 123 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB):1:1,000-1:5,000 Immunocytochemistry (IC): 1:100-1:1,000 Flow Cytometry (FCM): 1:2,000

Target information (Q14139)


Gene Name	UBE4A
About Gene Symbol	UBE4A (Ubiquitination Factor E4A) is a metabolic enzyme involved in ubiquitin-proteasome pathway. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.




Copyright ©2024 Genuin Biotechnologies LLC



Copyright ©2024 Genuin Biotechnologies LLC



Copyright ©2024 Genuin Biotechnologies LLC



Immunocytochemical staining of HepG2 cells with anti-ubiquitination factor E4 A antibody (Cat#ABG1717, 1:1,000). Nuclei were stained blue with DAPI; ubiquitination factor E4 A was stained magenta with Alexa Fluor® 647. Images were taken using Leica stellaris 5. Protein abundance based on laser Intensity and smart gain: Medium. Scale bar: 20 μ m.