

KD-Validated Anti-TACC3 Rabbit Monoclonal Antibody

Rabbit monoclonal antibody

ABG1720

Product Overview

Name	KD-Validated Anti-TACC3 Rabbit Monoclonal Antibody
Accession(Primary)	Q9Y6A5
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human TACC3
Species Reactivity	Human, Rat
Calculated MW	Predicted, 90 kDa , observed , 140 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB):1:1,000-1:5,000 Flow Cytometry (FCM):1:2,000

Target information (Q9Y6A5)

Gene Name	TACC3
About Gene Symbol	TACC3 (Transforming Acidic Coiled-Coil Containing Protein 3) is a formin actin nucleator involved in cytoskeletal organization. It is an important research target in cell biology, cancer metastasis, and cytoskeletal research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.

