

KD-Validated Anti-Tumor Protein P73 Rabbit Monoclonal Antibody

Rabbit monoclonal antibody


ABG1722

Product Overview


Name	KD-Validated Anti-Tumor Protein P73 Rabbit Monoclonal Antibody
Accession(Primary)	O15350
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human p73
Species Reactivity	Human
Calculated MW	Predicted, 70 kDa , observed , 70 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB):1:1,000-1:5,000 Flow Cytometry (FCM):1:2,000

Target information (O15350)


Gene Name	TP73
About Gene Symbol	TP73 (Tumor Protein P73) is a transcription factor involved in gene transcription and regulation. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Copyright ©2024 Genuin Biotechnologies LLC



Copyright ©2024 Genuin Biotechnologies LLC



Copyright ©2024 Genuin Biotechnologies LLC

