

KD-Validated Anti-TNFRSF12A Rabbit Monoclonal Antibody

Rabbit monoclonal antibody


ABG1723

Product Overview


Name	KD-Validated Anti-TNFRSF12A Rabbit Monoclonal Antibody
Accession(Primary)	Q9NP84
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human TWEAKR
Species Reactivity	Human
Calculated MW	Predicted, 14 kDa , observed , 18 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB):1:1,000-1:5,000 Flow Cytometry (FCM):1:2,000

Target information (Q9NP84)


Gene Name	TNFRSF12A
About Gene Symbol	TNFRSF12A (TNF Receptor Superfamily Member 12A) is a cytokine and immune signaling protein involved in adaptive immunity. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Copyright ©2024 Genuin Biotechnologies LLC



©2024 Genuin Biotechnol



Copyright ©2024 Genuin Biotechnologies LLC